


MEDITECH Long Term Care Ordering Guide
Respiratory Pathogen Testing - Influenza A & B/RSV/Covid-19/Magpix



**Note: Throat swabs are the least preferred as they are a suboptimal specimen type, and should only be used if the patient cannot tolerate a nasopharyngeal swab. Throat swabs must be collected with the “blue-top” Copan Universal Transport Medium kit.*

October to April (during Flu season): Influenza A & B/RSV and COVID-19 PCR testing will be done at an IH site for all specimens. “Magpix” testing will be performed at Kelowna General Hospital if specimen is PCR negative for COVID-19, Influenza A/B & RSV.

May to October (outside of Flu Season): Specimens will be tested for Covid-19. Influenza and RSV testing is not routinely available. “Magpix” testing will be performed at Kelowna General Hospital if specimen is PCR negative for COVID-19.

Note: If testing for Legionella is required, submit a lower respiratory specimen (sputum, endotracheal aspirate or bronchial wash) with a “Magpix” order. Legionella is not performed on nasopharyngeal samples or throat swabs.